

Mobile Shield

Bring a full 2.0mm Pb equivalency of radiation protection anywhere you need it with our x-ray protective mobile barriers

- Manufactured of the highest quality materials and craftsmanship
- 2.0mm Pb equivalency
- Transparent/clear upper section of x-ray protective Schott German-made shatterproof leaded glass
- Opaque lower section of solid steel with 2.0mm Pb Equiv. leaded x-ray protection.
- Frame color: Tan
- Locking caster system for easy mobility.
- Crated for secure shipment.
- Shipping worldwide by ground, air or sea freight
- Made in USA



Cat# PBMB-3630

Leaded Glass: 30" w x 36" h
Opaque Panel: 30" w x 39" h
Overall Size: 33" w x 75" h
Ships @ 450 lbs

